

# RRAGB Antibody



PACO03338

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

RRAGB

**Uniprot**

Q5VZM2

**Synonyms:**

RRAGB; Ras-related GTP-binding protein B; Rag B; RagB

**Immunogen:**

Synthesized peptide derived from the Internal region of human Rag B.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A