

Phospho-CEBPA (T226) Antibody



PACO03344

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

CEBPA

Uniprot

P49715

Synonyms:

CEBPA; CCAAT/enhancer-binding protein alpha; C/EBP alpha

Immunogen:

Synthesized peptide derived from human C/EBP alpha around the phosphorylation site of T226.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A