

TRIM16 Antibody



PACO03359

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

TRIM16

Uniprot

O95361

Synonyms:

TRIM16; EBBP; Tripartite motif-containing protein 16; Estrogen-responsive B box protein

Immunogen:

Synthesized peptide derived from the N-terminal region of human TRIM16.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A