

# SIX2 Antibody



PACO03366

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

SIX2

**Uniprot**

Q9NPC8

**Synonyms:**

SIX2; Homeobox protein SIX2; Sine oculis homeobox homolog 2

**Immunogen:**

Synthesized peptide derived from the Internal region of human Six2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A