

BLCAP Antibody



PACO03370

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

BLCAP

Uniprot

P62952

Synonyms:

BLCAP; BC10; Bladder cancer-associated protein; Bladder cancer 10 kDa protein; Bc10

Immunogen:

Synthesized peptide derived from the Internal region of human Bc10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A