

# ALPK2 Antibody



PACO03401

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ALPK2

**Uniprot**

Q86TB3

**Synonyms:**

ALPK2; HAK; Alpha-protein kinase 2; Heart alpha-protein kinase

**Immunogen:**

Synthesized peptide derived from the Internal region of human Alpha-protein Kinase 2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A