

DHRS1 Antibody



PACO03402

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

DHRS1

Uniprot

Q96LJ7

Synonyms:

DHRS1; Dehydrogenase/reductase SDR family member 1

Immunogen:

Synthesized peptide derived from the Internal region of human DHRS1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A