

RPL39L Antibody



PACO03406

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

RPL39L

Uniprot

Q96EH5

Synonyms:

RPL39L; RPL39L1; 60S ribosomal protein L39-like; 60S ribosomal protein L39-2

Immunogen:

Synthesized peptide derived from the Internal region of human Ribosomal Protein L39L.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A