

# Phospho-ADD1 (T445) Antibody



PACO03410

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ADD1

**Uniprot**

P35611

**Synonyms:**

ADD1; ADDA; Alpha-adducin; Erythrocyte adducin subunit alpha

**Immunogen:**

Synthesized peptide derived from human Adducin alpha around the phosphorylation site of T445.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A