

# OR4B1 Antibody



PACO03414

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR4B1

**Uniprot**

Q8NGF8

**Synonyms:**

OR4B1; Olfactory receptor 4B1; OST208; Olfactory receptor OR11-106

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4B1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A