

OR10J5 Antibody



PACO03431

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR10J5

Uniprot

Q6IFH2

Synonyms:

Olfactory receptor 10J5

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 10J5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A