

PODN Antibody



PACO03432

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

PODN

Uniprot

Q7Z5L7

Synonyms:

PODN; SLRR5A; Podocan

Immunogen:

Synthesized peptide derived from the C-terminal region of human Podocan.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A