

COL4A3 Antibody



PACO03435

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

COL4A3

Uniprot

Q01955

Synonyms:

COL4A3; Collagen alpha-3(IV) chain; Goodpasture antigen

Immunogen:

Synthesized peptide derived from the Internal region of human COL4A3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A