

# CSTL1 Antibody



PACO03439

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

CSTL1

**Uniprot**

Q9H114

**Synonyms:**

CSTL1; Cystatin-like 1; RCET11

**Immunogen:**

Synthesized peptide derived from the Internal region of human Cystatin L1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A