

# COL7A1 Antibody



PACO03445

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

COL7A1

**Uniprot**

Q02388

**Synonyms:**

COL7A1; Collagen alpha-1(VII) chain; Long-chain collagen; LC collagen

**Immunogen:**

Synthesized peptide derived from the Internal region of human COL7A1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A