

# COL12A1 Antibody



PACO03447

---

## Product Information

<b>Size:</b> 50ug	<b>Protein Background:</b>
<b>Reactivity:</b> Human	<b>Gene ID:</b> COL12A1
<b>Source:</b> Rabbit	<b>Uniprot</b> Q99715
<b>Isotype:</b> IgG	<b>Synonyms:</b> COL12A1; COL12A1L; Collagen alpha-1(XII) chain
<b>Applications:</b> ELISA, IHC, IF	<b>Immunogen:</b> Synthesized peptide derived from the Internal region of human COL12A1.
<b>Recommended dilutions:</b> IHC:1:100-1:300, IF:1:200-1:1000	<b>Storage:</b> Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A