

# COL19A1 Antibody



PACO03453

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

COL19A1

**Uniprot**

Q14993

**Synonyms:**

COL19A1; Collagen alpha-1(XIX) chain; Collagen alpha-1(Y) chain

**Immunogen:**

Synthesized peptide derived from the Internal region of human COL19A1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A