

GLYCTK Antibody



PACO03457

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

GLYCTK

Uniprot

Q8IVS8

Synonyms:

GLYCTK; HBEBP4; LP5910; Glycerate kinase; HBeAg-binding protein 4

Immunogen:

Synthesized peptide derived from the Internal region of human Glycerate Kinase.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A