

# OR2Y1 Antibody



PACO03464

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR2Y1

**Uniprot**

Q8NGV0

**Synonyms:**

OR2Y1; Olfactory receptor 2Y1; Olfactory receptor OR5-2

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 2Y1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A