

# OR13C8 Antibody



PACO03495

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

OR13C8

**Uniprot**

Q8NGS7

**Synonyms:**

OR13C8; Olfactory receptor 13C8

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 13C8.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A