

OR13C3 Antibody



PACO03496

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR13C3

Uniprot

Q8NGS6

Synonyms:

OR13C3; Olfactory receptor 13C3; Olfactory receptor OR9-8

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 13C3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A