

PAGE3 Antibody



PACO03504

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

PAGE3

Uniprot

Q5JUK9

Synonyms:

PAGE3; GAGED1; Putative G antigen family D member 1; Prostate-associated gene 3 protein; PAGE-3

Immunogen:

Synthesized peptide derived from the Internal region of human PAGE-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A