

TMEM37 Antibody



PACO03512

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

TMEM37

Uniprot

Q8WXS4

Synonyms:

TMEM37; PR; Voltage-dependent calcium channel gamma-like subunit; Neuronal voltage-gated calcium channel gamma-like subunit; Transmembrane protein 37

Immunogen:

Synthesized peptide derived from the N-terminal region of human TMEM37.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A