

CST2 Antibody



PACO03531

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

CST2

Uniprot

P09228

Synonyms:

CST2; Cystatin-SA; Cystatin-2; Cystatin-S5

Immunogen:

Synthesized peptide derived from the Internal region of human Cystatin SA.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A