

CST6 Antibody



PACO03532

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

CST6

Uniprot

Q15828

Synonyms:

CST6; Cystatin-M; Cystatin-6; Cystatin-E

Immunogen:

Synthesized peptide derived from the Internal region of human Cystatin E/M.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A