

# CNTN4 Antibody



PACO03545

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:****Gene ID:**

CNTN4

**Uniprot**

Q8IWW2

**Synonyms:**

CNTN4; Contactin-4; Brain-derived immunoglobulin superfamily protein 2; BIG-2

**Immunogen:**

Synthesized peptide derived from the Internal region of human Contactin 4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A