

CYP2A7 Antibody



PACO03553

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

CYP2A7

Uniprot

P20853

Synonyms:

CYP2A7; Cytochrome P450 2A7; CYP11A7; Cytochrome P450 IIA4

Immunogen:

Synthesized peptide derived from the Internal region of human CYP2A6/7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A