

# DOCK1 Antibody



PACO03576

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

DOCK1

**Uniprot**

Q14185

**Synonyms:**

DOCK1; Deducator of cytokinesis protein 1; 180 kDa protein downstream of CRK; DOCK180

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human DOCK 180.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A