

DOCK2 Antibody



PACO03577

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

DOCK2

Uniprot

Q92608

Synonyms:

DOCK2; KIAA0209; Deducator of cytokinesis protein 2

Immunogen:

Synthesized peptide derived from the Internal region of human DOCK 2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A