

EGR3 Antibody



PACO03591

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

EGR3

Uniprot

Q06889

Synonyms:

EGR3; PILOT; Early growth response protein 3; EGR-3; Zinc finger protein pilot

Immunogen:

Synthesized peptide derived from the C-terminal region of human Egr-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A