

# CYP4Z1 Antibody



PACO03605

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

CYP4Z1

**Uniprot**

Q86W10

**Synonyms:**

CYP4Z1; Cytochrome P450 4Z1; CYP4Z1

**Immunogen:**

Synthesized peptide derived from the Internal region of human CYP4Z1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A