

OR8U1/OR8U8/OR8U9 Antibody



PACO03654

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR8U1/OR8U8/OR8U9

Uniprot

Q8NH10/P0C7N1/P0C7N5

Synonyms:

OR8U1; Olfactory receptor 8U1; OR8U8; Olfactory receptor 8U8; OR8U9; Olfactory receptor 8U9

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 8U1/8/9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A