

GPR116 Antibody



PACO03659

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

ADGRF5

Uniprot

Q8IZF2

Synonyms:

GPR116; KIAA0758; Probable G-protein coupled receptor 116

Immunogen:

Synthesized peptide derived from the C-terminal region of human GPR116.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A