

LPHN1 Antibody



PACO03669

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

ADGRL1

Uniprot

O94910

Synonyms:

LPHN1; KIAA0821; LEC2; Latrophilin-1; Calcium-independent alpha-latrotoxin receptor 1; CIRL-1; Lectomedin-2

Immunogen:

Synthesized peptide derived from the Internal region of human Latrophilin-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A