

MYH15 Antibody



PACO03674

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

MYH15

Uniprot

Q9Y2K3

Synonyms:

MYH15; KIAA1000; Myosin-15; Myosin heavy chain 15

Immunogen:

Synthesized peptide derived from the Internal region of human MYH15.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A