

# FOXI1 Antibody



PACO03675

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

FOXI1

**Uniprot**

Q12951

**Synonyms:**

FOXI1; FKHL10; FREAC6; Forkhead box protein I1; Forkhead-related protein FKHL10; Forkhead-related transcription factor 6; FREAC-6; Hepatocyte nuclear factor 3 forkhead homolog 3; HFH-3; HNF-3/fork-head homolog 3

**Immunogen:**

Synthesized peptide derived from the Internal region of human FoxI1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A