

ZC3H4 Antibody



PACO03689

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

ZC3H4

Uniprot

Q9UPT8

Synonyms:

ZC3H4; C19orf7; KIAA1064; Zinc finger CCCH domain-containing protein 4

Immunogen:

Synthesized peptide derived from the Internal region of human ZC3H4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A