

KCNH3 Antibody



PACO03700

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

KCNH3

Uniprot

Q9ULD8

Synonyms:

KCNH3; KIAA1282; Potassium voltage-gated channel subfamily H member 3; Brain-specific eag-like channel 1; BEC1; Ether-a-go-go-like potassium channel 2; ELK channel 2; ELK2; Voltage-gated potassium channel subunit Kv12.2

Immunogen:

Synthesized peptide derived from the Internal region of human KCNH3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A