

RASD2 Antibody



PACO03716

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

RASD2

Uniprot

Q96D21

Synonyms:

RASD2; TEM2; GTP-binding protein Rhes; Ras homolog enriched in striatum; Tumor endothelial marker 2

Immunogen:

Synthesized peptide derived from the C-terminal region of human Dexas2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A