

LDOC1 Antibody



PACO03733

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

LDOC1

Uniprot

O95751

Synonyms:

LDOC1; BCUR1; Protein LDOC1; Leucine zipper protein down-regulated in cancer cells

Immunogen:

Synthesized peptide derived from the Internal region of human LDOC1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A