

ZKSCAN5 Antibody



PACO03735

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

ZKSCAN5

Uniprot

Q9Y2L8

Synonyms:

ZKSCAN5; KIAA1015; ZFP95; Zinc finger protein with KRAB and SCAN domains 5; Zinc finger protein 95 homolog; Zfp-95

Immunogen:

Synthesized peptide derived from the Internal region of human ZKSCAN5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A