

RAB38 Antibody



PACO03736

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

RAB38

Uniprot

P57729

Synonyms:

RAB38; Ras-related protein Rab-38; Melanoma antigen NY-MEL-1

Immunogen:

Synthesized peptide derived from the Internal region of human Rab 38.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A