

# Phospho-ALOX5 (S272) Antibody



PACO03747

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

ALOX5

**Uniprot**

P09917

**Synonyms:**

ALOX5; LOG5; Arachidonate 5-lipoxygenase; 5-LO; 5-lipoxygenase

**Immunogen:**

Synthesized peptide derived from human 5-LO around the phosphorylation site of S272.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A