

ACKR1 Antibody



PACO03768

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

ACKR1

Uniprot

Q16570

Synonyms:

DARC; FY; GPD; Duffy antigen/chemokine receptor; Fy glycoprotein; GpFy; Glycoprotein D; Plasmodium vivax receptor; CD antigen CD234

Immunogen:

Synthesized peptide derived from the N-terminal region of human Duffy.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A