

BRWD3 Antibody



PACO03772

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

BRWD3

Uniprot

Q6RI45

Synonyms:

BRWD3; Bromodomain and WD repeat-containing protein 3

Immunogen:

Synthesized peptide derived from the C-terminal region of human BRWD3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A