

# PART1 Antibody



PACO03791

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

PART1

**Uniprot**

Q9NPD0

**Synonyms:**

Prostate-specific and androgen-regulated protein

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human PART-1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A