

GALR1 Antibody



PACO03793

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:**Gene ID:**

GALR1

Uniprot

P47211

Synonyms:

GALR1; GALNR; GALNR1; Galanin receptor type 1; GAL1-R; GALR-1

Immunogen:

Synthesized peptide derived from the Internal region of human GALR1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A