

TAS2R43 Antibody



PACO03799

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

TAS2R43

Uniprot

P59537

Synonyms:

TAS2R43; Taste receptor type 2 member 43; T2R43; Taste receptor type 2 member 52; T2R52

Immunogen:

Synthesized peptide derived from the Internal region of human T2R43.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A