

LRP10 Antibody



PACO03808

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

LRP10

Uniprot

Q7Z4F1

Synonyms:

LRP10; MSTP087; SP220; Low-density lipoprotein receptor-related protein 10; LRP-10

Immunogen:

Synthesized peptide derived from the Internal region of human LRP10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A