

Phospho-GATA2 (S401) Antibody



PACO03826

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

GATA2

Uniprot

P23769

Synonyms:

GATA2; Endothelial transcription factor GATA-2; GATA-binding protein 2

Immunogen:

Synthesized peptide derived from human GATA-2 around the phosphorylation site of S401.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A